
**Seprafilm CE insert**

Proof 3 r92 02.12.20 op:rs area: box language(s): multiple  
code: GEN job: 510677 10%\* x 24\* 4564 (02/20)

Insert is 100% of actual size

**BARRIÈRE ANTI-ADHÉRENCES**  
**ADHÉSIOBARRIÈRE**  
**ADHESIEBARRIÈRE**  
**ADHERANSEBARRIÈRE**  
**BARRIERA ANTIADHERENZE**  
**BARRERA ANTIADHERENCIA**  
**BARREIRA ANTIADHERÊNCIA**  
**ΦΑΓΜΑ ΣΥΜΦΥΣΣΕΩΝ**

Seprafilm and design is a registered trademark of Genzyme Corporation, Cambridge, MA, USA.

Genzyme is a registered trademark of Genzyme Corporation, Cambridge, MA, USA.

Genzyme Europe B.V.  
Gooimeer 10  
1411 DD Naarden  
The Netherlands  
Tel: +31-35-699-1200  
Fax: +31-35-694-3214

In patients undergoing surgery for ovarian, primary peritoneal or fallopian tube malignancies, Seprafilm use has been reported as having an increased risk of intrabdominal fluid collection and/or abscess, particularly when extensive debulking surgery was required.

- CAUTIONS**
- Do not use if foil pouch is damaged or appears to have been tampered with or opened prior to time of use.
  - The safety and effectiveness of Seprafilm in combination with other adhesion prevention products and/or in other surgical procedures not with the abdominopelvic cavity have not been established in clinical studies.
  - No controlled clinical studies have been conducted in patients with active infections.
  - Foreign body reactions may occur, as with any implanted material, most surgical adjuncts, but have been rarely reported during clinical use.
  - Seprafilm has limited controlled clinical trial information in patients with abdominopelvic malignancy.
  - No pre-clinical reproductive studies have been conducted. No clinical studies have been conducted in women who become pregnant in the first month after application of Seprafilm. Therefore, avoiding pregnancy during the first complete menstrual cycle after the use of Seprafilm should be considered.
- ADVERSE REACTIONS**
- During post-marketing experience, the following adverse events have been reported in patients receiving Seprafilm Adhesion Barrier: abscess, anastomotic leak, fistula, foreign body (fibrotic) reaction, hypersensitivity, inflammation, infectious fluid collection, peritonitis, post-procedural wound infection/wound dehiscence, sepsis.

CE 2797  
4564 (02/20)

- PMS 534
- PMS 356
- Black

**DESCRIPTION**

Seprafilm® Adhesion Barrier (Seprafilm) is a sterile, bioresorbable, translucent membrane composed of two chemically modified anionic polysaccharides, sodium hyaluronate (HA) and carboxymethylcellulose (CMC).

**INDICATION**

Seprafilm is intended as an adjunct in abdominal and pelvic surgery for reducing the incidence, extent and severity of postoperative adhesions at the site of placement, and to reduce adhesive small bowel obstruction when placed in the abdomen.

**INTENDED USE**

Seprafilm should be applied to sites of potentially adhesogenic tissue and organ structures in the abdominopelvic cavity to serve as a temporary barrier separating opposing tissue surfaces.

**CONTRAINDICATIONS**

Seprafilm is contraindicated in patients with a history of hypersensitivity to Seprafilm and/or to any component of Seprafilm.

**WARNINGS**

Seprafilm must be used according to the instructions for use. Read instructions prior to use. Seprafilm is supplied sterile and should not be resterilized. The membrane is for single use only. Every opened and unused Seprafilm pouch must be discarded.

Seprafilm is not recommended to be wrapped directly around a fresh anastomotic suture or staple line of the intestine. Clinical trial data on Seprafilm indicate that such use may result in an increased risk of anastomotic leak-related events (fistula, abscess, leak, sepsis and peritonitis). The incidence of these events was not affected when Seprafilm was placed elsewhere in the abdomen.

In patients undergoing surgery for ovarian, primary peritoneal or fallopian tube malignancies, Seprafilm use has been reported as having an increased risk of intrabdominal fluid collection and/or abscess, particularly when extensive debulking surgery was required.

- Do not use if foil pouch is damaged or appears to have been tampered with or opened prior to time of use.
- The safety and effectiveness of Seprafilm in combination with other adhesion prevention products and/or in other surgical procedures not with the abdominopelvic cavity have not been established in clinical studies.
- No controlled clinical studies have been conducted in patients with active infections.
- Foreign body reactions may occur, as with any implanted material, most surgical adjuncts, but have been rarely reported during clinical use.
- Seprafilm has limited controlled clinical trial information in patients with abdominopelvic malignancy.
- No pre-clinical reproductive studies have been conducted. No clinical studies have been conducted in women who become pregnant in the first month after application of Seprafilm. Therefore, avoiding pregnancy during the first complete menstrual cycle after the use of Seprafilm should be considered.

**CAUTIONS**

- Do not use if foil pouch is damaged or appears to have been tampered with or opened prior to time of use.
- The safety and effectiveness of Seprafilm in combination with other adhesion prevention products and/or in other surgical procedures not with the abdominopelvic cavity have not been established in clinical studies.
- No controlled clinical studies have been conducted in patients with active infections.
- Foreign body reactions may occur, as with any implanted material, most surgical adjuncts, but have been rarely reported during clinical use.
- Seprafilm has limited controlled clinical trial information in patients with abdominopelvic malignancy.
- No pre-clinical reproductive studies have been conducted. No clinical studies have been conducted in women who become pregnant in the first month after application of Seprafilm. Therefore, avoiding pregnancy during the first complete menstrual cycle after the use of Seprafilm should be considered.

**ADVERSE REACTIONS**

During post-marketing experience, the following adverse events have been reported in patients receiving Seprafilm Adhesion Barrier: abscess, anastomotic leak, fistula, foreign body (fibrotic) reaction, hypersensitivity, inflammation, infectious fluid collection, peritonitis, post-procedural wound infection/wound dehiscence, sepsis.

**STORAGE INSTRUCTIONS**

Seprafilm should be stored at 2-30°C.

**HOW SUPPLIED**

Seprafilm is packed in a Tyvek® holder under a plastic sleeve and packed in an outer, sealed, foil pouch. The contents of the foil pouch are sterilized by gamma radiation of 25-45 kGy.

Refer to package label for film size and quantity.

**DIRECTIONS FOR USE**

- Open the foil pouch immediately prior to application and drop the inferior sterile polyolefin sleeve containing Seprafilm Adhesion Barrier on the dry sterile field.
- Carefully remove the Tyvek® holder from the plastic sleeve.
- Keep membrane dry in the holder prior to application.
- While membrane remains in the holder, cut to desired size with scissors.
- Expose 2 cm of membrane through open end of the holder.

**Placement**

Apply Seprafilm immediately prior to closure of the abdominopelvic cavity:

- Ensure surgical field is dry and aspirate excess fluid.
- Handle membrane carefully with dry instruments and/or gloves.
- Avoid contact with tissue surfaces until directly at site of application. If contact occurs, moderate application of a standard irrigating solution may be used to gently dislodge membrane from unintended tissue surfaces.
- When necessary, facilitate entry in abdominal and pelvic cavity by slightly curving the holder.
- Extend placement up to 8 cm beyond margins of direct trauma to achieve adequate coverage with the membrane.
- Let exposed membrane adhere to desired position on the tissue or organ by gently pressing down with a dry glove or instrument, while withdrawing the holder.
- When necessary, moisten membrane lightly with irrigating solution (1-2 ml) to facilitate moulding of membrane along the tissue or around the organ contours.
- Up to 10 (13 cm x 15 cm) membranes per patient have been used in controlled clinical studies for the abdominopelvic indication. Allow overlap of 2-3 cm between individual membranes to ensure coverage.

**CONTRAINDICATIONS**

Seprafilm is contraindicated in patients with a history of hypersensitivity to Seprafilm and/or to any component of Seprafilm.

In patients undergoing surgery for ovarian, primary peritoneal or fallopian tube malignancies, Seprafilm use has been reported as having an increased risk of intrabdominal fluid collection and/or abscess, particularly when extensive debulking surgery was required.

- Do not use if foil pouch is damaged or appears to have been tampered with or opened prior to time of use.
- The safety and effectiveness of Seprafilm in combination with other adhesion prevention products and/or in other surgical procedures not with the abdominopelvic cavity have not been established in clinical studies.
- No controlled clinical studies have been conducted in patients with active infections.
- Foreign body reactions may occur, as with any implanted material, most surgical adjuncts, but have been rarely reported during clinical use.
- Seprafilm has limited controlled clinical trial information in patients with abdominopelvic malignancy.
- No pre-clinical reproductive studies have been conducted. No clinical studies have been conducted in women who become pregnant in the first month after application of Seprafilm. Therefore, avoiding pregnancy during the first complete menstrual cycle after the use of Seprafilm should be considered.

**CAUTIONS**

- Do not use if foil pouch is damaged or appears to have been tampered with or opened prior to time of use.
- The safety and effectiveness of Seprafilm in combination with other adhesion prevention products and/or in other surgical procedures not with the abdominopelvic cavity have not been established in clinical studies.
- No controlled clinical studies have been conducted in patients with active infections.
- Foreign body reactions may occur, as with any implanted material, most surgical adjuncts, but have been rarely reported during clinical use.
- Seprafilm has limited controlled clinical trial information in patients with abdominopelvic malignancy.
- No pre-clinical reproductive studies have been conducted. No clinical studies have been conducted in women who become pregnant in the first month after application of Seprafilm. Therefore, avoiding pregnancy during the first complete menstrual cycle after the use of Seprafilm should be considered.

**ADVERSE REACTIONS**

During post-marketing experience, the following adverse events have been reported in patients receiving Seprafilm Adhesion Barrier: abscess, anastomotic leak, fistula, foreign body (fibrotic) reaction, hypersensitivity, inflammation, infectious fluid collection, peritonitis, post-procedural wound infection/wound dehiscence, sepsis.

**All others contact:**  
Genzyme Biosurgery  
4620 Kingsgate  
Cascade Way  
Oxford Business Park South  
Oxford  
OX4 2SU  
Tel: +44 (0)1865 405 200  
Fax: +44 (0)1865 774 172

\*Tyvek® is a registered trademark of the DuPont Company, Wilmington, DE, U.S.A.

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